


NTEGRA THz

THz near-field imaging system (sSNOM).

THz near-field microscope has been developed together with a new generation of NTEGRA system. We utilize the high-power laser to generate a stable THz pulse train. Single-cycle THz pulses with THz electric fields and a 0.15 – 5 THz spectral width can be generated with the photoconductive antenna system. The THz beam is focused down to a $\sim 600 \mu\text{m}$ spot under the AFM tip which yield scattering-type THz near-field signal with $< 50 \text{ nm}$ spatial resolution. The THz near-field system is also equipped with built-in heating stage and purge box, capable of studying high-temperature THz near-field spectroscopy with different gaseous environment.



Applications

Specifications

SCHEMATIC OF THE NEAR FIELD OPTICAL-PUMP THZ-PROBE (N-OPTP) SETUP




Figure 1: Schematic of the near field optical-pump THz-probe (n-OPTP) setup, equally capable of performing optical pump near-field THz emission (n-OPTE) experiments. ①: 800 nm, 300 mW pump pulses pass through the ITO and collimate with the THz beam ②. Optical pump ① and THz probe ② are focused onto the AFM tip apex using an off-axis parabolic mirror. ③: THz gate (detection) beam. Tip scattered THz signals can be mapped out in the time domain by changing the time delay between ② and ③ (t_1). Changing the delay between ① and ③ (t_2) probes the photo-excited ultrafast dynamics of the sample. When t_1 is fixed at the peak position of the scattered THz electric field while t_2 is varied, n-OPTE (with THz probe ② blocked) and n-OPTP (with THz probe ② unblocked) can be performed.

SPECIFICATIONS

Far-field THz performance:

Excitation Laser	Laser central wavelength 800 nm +/- 50 nm, pulse length < 120 fs, average power > 100 mW; time synchronized with the THz pulses
Measurement Modalities	Transmission & Reflection
Average Optical Power on Transmitter	> 15 mW
Average Optical Power on Receiver	> 15 mW
Bias Voltage on Transmitter	± 60V square wave
Terahertz Peak Measured Photocurrent	> 2 nA high BW; > 30 nA high signal
Terahertz Spectrum Bandwidth	> 5 THz high BW; > 4 THz high signal
Power Spectrum Dynamic Range	> 70 dB high BW; > 80 dB high signal
Typical Scan Time	2-5 min

Dry air purge box

Imaging raster scan hardware and software

Lock-in amplifier and Low-Noise current Amplifier

Far-field THz performance:

The basic laser and scanning scheme are the same as the far-field setup.

Near-field THz spectral range:	0.3- 2 THz
Near-field THz signal to noise:	S2 20:1 at 30ms lock-in integration time
Typical Scan Time	15-60 min

Dry air purge box

Near-field Image taken in integrated NT-MDT AFM software

Optical pump THz probe compatible

THz near-field emission microscopy compatible

APPLICATION EXAMPLE

Using the THz near-field system, nano-imaging of micrometer size graphene samples on SiO₂ (300 nm)/Si substrate were performed. The graphene was mechanically exfoliated onto PDMS (Polydimethylsiloxane) and dry transferred onto a prepatterned Au lead (30 nm thick). Here the Au lead is used for electrostatic gating and serves as a reference for graphene imaging. Figure 2a,b show the simultaneously collected AFM topography and near-field images with a size of 10 × 10 μm, respectively. As shown in Figure 2b, there is clearly a large optical contrast between single layer graphene (SLG, marked as 1) and the SiO₂ substrate (marked as 0). However, no obvious contrasts between SLG and multilayer graphene can be observed (the multilayers are marked by the number of layers as 2, 3, 4, and 5+). Moreover, a comparison between single-layer graphene and a 30 nm gold film revealed comparable THz near-field signals (see Figure 2c,d). These observations suggest that monolayer graphene has close-to- perfect near-field THz reflection.³¹ This near-unity THz near- field reflectivity is found to be independent of the carrier density when the back-gate voltage is swept between ±30 V at room temperature (RT). To compare the THz near-field spectra between SiO₂ and graphene, we plot the time domain signal in Figure 2e and the Fourier-transformed THz spectrum in Figure 2f. The normalized THz nanospectrum $S_2(\text{graphene})/S_2(\text{SiO}_2)$ reveals that the near-field signal of graphene is 3–4× larger than that on the bare SiO₂ over the 0.2–1 THz frequency range. In the inset of Figure 2f, the unnormalized spectra of graphene and SiO₂ are plotted. A dip at about 0.75 THz can be observed, which was reported in previous THz near-field experiments. Since this dip is absent in the far-field THz measurements using the same optics but excluding the involvement of the AFM tip, it is considered to be caused by the complex tip-light interactions and can be an interesting point for further investigation.




Figure 2: (a, b) AFM topography and THz near-field (S_2) mapping of graphene on SiO₂, respectively. The numbers of graphene layers are marked in (a), with bare SiO₂ marked as 0. (c, d) AFM and THz near-field (S_2) images of a SLG with a gold electrode. The near-field signal in graphene is comparable to that on the thin gold films. (e) Near field THz-TDS signal of SiO₂ (black) and graphene (red). (f) Normalized graphene THz near- field spectrum (to SiO₂). The inset shows a Fast Fourier Transform (FFT) of (e), which is the unnormalized S_2 spectra of graphene (red) and SiO₂ (black).

For more details please refer to Zhang, Jiawei, Chen, Xinzhong, et al. "Terahertz Nano-imaging of graphene." *ACS Photonics* (2018).



+31 641 494682

info@ntmdt.nl

Hoenderparkweg 96B, 7335 GX
Apeldoorn, the Netherlands

Products

NTEGRA Prima

NTEGRA MFM

NTEGRA Aura

NTEGRA SNOM

NTEGRA Spectra

NTEGRA THz

Company

ABOUT US

EVENTS

PARTNERS

CONTACTS